



**SPECIALTY MEDICAL PRODUCTS, INC.**

**Billing:** 50 Pennsylvania Ave Malvern, PA 19382

**Shipping:** 187 Pennsylvania Ave Malvern, PA 19382

**Phone:** (888) 698-1522      **Fax (610) 644-3992**

**Email:** [info@smpcares.com](mailto:info@smpcares.com)

**Website:** [www.smpcares.com](http://www.smpcares.com)



## The Merry Walker® Ambulation Device: INSTITUTIONAL

The Merry Walker® Ambulation Device is a walker/chair combination, allowing a person who would normally be placed in a wheelchair to be placed into a Merry Walker, allowing independent and safe walking. End users will be able to ambulate independently and safely.

The Merry Walker is fabricated from one inch, 16 gauge steel, is bottom weighted to prevent unwanted tipping, and the unit folds easily for storage and transport.

Accessories and replacement parts are also available (see below). The Merry Walker® Ambulation Device comes in three models, and the INSTITUTIONAL model is featured below:

Merry Walker® Institutional - Used in institutions such as hospitals and nursing homes. Base 27" wide. Weight capacity up to 300 lb.

### Benefits:

- Standard pick-up walker
- Restraint-free or justified/enabler
- Returns independence to people
- Enhances self-rehabilitation
- Reduces the need to use two people to walk one person.
- Keeps people ambulatory and allows them to walk independently, not relying on assistance.
- Increases appetites.
- Decreases learned dependency.
- Encourages appropriate sleep patterns.
- Addresses the dignity issue as people are up and walking independently.
- Safety strap for safety/optional Vee Belt for Alzheimer's Disease and related dementias diagnosed people to use.

### Features:

- Black powder coated tubular steel/washable.
- Sized by height and weight of each person.
- Cushioned, washable seat. (Gel Seat sold separately)
- Side pouch for personal items.
- All fabrics washable.
- Swell latch lock, adjustable tension.