

Invia[®] Liberty[™] NPWT System


THERAPY CONFIDENCE MADE SIMPLE

Therapy Confidence

- **Reliable** negative pressure delivery and control at the wound site¹
- **Intelligent** blockage detection¹ from pump to dressing
- **Secure** canister and tubing connection² with an audible click

Simplified Management


- **Intuitive** pump operation and troubleshooting instructions
- **Compact** design for patient comfort and mobility
- **Clinical flexibility** with comprehensive dressings portfolio



Medela provides efficient solutions to support clinicians with Negative Pressure Wound Therapy management in all patient care settings.

Invia Liberty delivers therapy **confidence by actively controlling set pressure at the wound site**¹. The pump quickly adapts the air flow cycle to the amount of fluid in the tubing for optimal system performance.

Click-in canister attachment and **one-touch canister release** button allow for easy canister exchange. **Secure tubing attachment** is guaranteed with the new Invia Quick-connector feature².








Intuitive and quiet pump operation promotes full integration in the patient's daily life, while **patient comfort and mobility** are facilitated by Invia Liberty's compact design and convenient carrying case.

Key features at a glance

- Active exudate management and pressure control from the pump to the wound site
- Pressure options from -40 to -200 mmHg including standard -125 mmHg
- On-screen air leak indicator
- New custom air leak modes for normal and large wounds
- Easy and secure tubing attachment with the Quick-connector
- One-touch canister release button
- Choice of 300 ml and 800 ml disposable canisters with solidifier



Technical data

 - 40 mmHg to -200 mmHg	 > 10 Hours	Modes: Constant and Intermittent Tubing: Double Lumen Tubing, Y-Connector Canister: 300 and 800 ml
 Low flow 5 l/min.	 1 kg 2.2 lbs	 H x W x D 95 x 170 x 150 mm 3.74 x 6.69 x 5.91 inch
		2 year warranty

Find out more about the new Invia Liberty system at www.medela-healthcare.us

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

¹ Performance data on file at Medela AG
² Usability Validation data on file at Medela AG

Medical Vacuum Technology
for Healthcare Professionals

Please contact us or your local Medela representative for details.

Local contact:



SPECIALTY MEDICAL PRODUCTS, INC.

Billing: 50 Pennsylvania Ave Malvern, PA 19382
Shipping: 187 Pennsylvania Ave Malvern, PA 19382
Phone: (888) 698-1522 Fax (610) 644-3992
Email: info@smpcares.com
Website: www.smpcares.com

Medela and Invia are registered in the US Patent and Trademark Office. Liberty is a trademark of Medela Holding AG, Switzerland. No trademark registrations in Canada.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

USA
Medela LLC
1101 Corporate Drive
McHenry, Illinois 60050
USA
Phone +1 877 735 1626
Fax +1 815 307 8942
info-healthcare@medela.com
www.medela-healthcare.us

Canada
Medela Canada Inc.
4160 Sladeview Cres., Unit #8
Mississauga, Ontario, L5L 0A1
Canada
Phone +1 800 435 8316
Fax +1 800 995 7867
info@medela.ca
www.medela.ca