




SPECIALTY MEDICAL PRODUCTS, INC.
50 PENNSYLVANIA AVE; MALVERN, PA 19355
(O) 610-644-1370 (FX) 610-644-3992

VasoPress



SPECIFICATIONS:

- Weight:** 6.4 pounds
- Pressure Range:** 40 mmHg (Leg); 80 mmHg (Foot)
- Power:** 200 ma, 15W
- Time:** Inflation: 12 seconds (approximate)
Deflation: 48 seconds (approximate)
- Indicators:** Power ON, Low Pressure and Continuous Pressure
- Alarms:** Audible, intermittent beep and visible LED light
- Voltage:** 120V ± 10%, 60 Hz, - Also Available in 220V
- Safety:** UL, CE, TUV

ORDERING INFORMATION:

- VP 500:** Vaso Press Pump
- VP 501P:** Vaso Press Calf Garment, Pediatric, for calf sizes up to 14 inches
- VP 501M:** Vaso Press Calf Garment, Medium, for calf sizes up to 18 inches
- VP 501L:** Vaso Press Calf Garment, Large, for calf sizes up to 24 inches
- VP 501B:** Vaso Press Calf Garment, Bariatric, for calf sizes up to 32 inches
- VP 520:** Vaso Press Foot Garment, Fits Right or Left Foot
- VP 530M:** Vaso Press Thigh Garment, Medium, for thigh sizes up to 29 inches
- VP 530L:** Vaso Press Thigh Garment, Large, for thigh sizes up to 36 inches
- VP 530B:** Vaso Press Thigh Garment, Bariatric, for thigh sizes up to 42 inches
- All garments are packaged as pairs, 10 pairs/case
- VP 501I:** Vaso Press Insert Piece, 5 pairs/case
- VP 510:** Vaso Press Tube Set, 60 inches, 5 pairs/case
- VP 510L:** Vaso Press Tube Set, 120 inches long, 5 pairs/case

- Flam E et al: Femoral Vein Blood Velocities With Intermittent Compression Systems: Implications For DVT Management. Poster Presentation, Feb 24, 1992, American Academy of Orthopedic Surgeons, Washington, D.C.
- Flam E et al: Blood-Flow Augmentation of Intermittent Pneumatic Compression Systems Used for the Prevention of Deep Vein Thrombosis Prior to Surgery. The American Journal of Surgery, March 1996.
- Salzman E.W. et al: Effect of Optimization of Hemodynamics on Fibrinolytic Activity and Antithrombotic Efficacy of External Pneumatic Calf Compression; Annals of Surgery, 1987; 206, No. 5:636-641.
- Healthcare Product Comparison System, ECRI, December 2001.
- Clyde M, Nelson B, Lisek, E. et al: The Comparison of Three Pneumatic Compression Devices for the Prevention of Deep Vein Thrombosis in Total Joint Replacement Patients. Poster Presentation, Feb 1995, American Academy of Orthopedic Surgeons.
- Camerota AJ, Katz ML, White JV: Why does prophylaxis with external pneumatic compression for deep vein thrombosis fail? Am J Surg 164:3;265-268; September 1992.
- Fedullo P, Brewer N, et al: Comparison of Thigh and Calf Length Intermittent Pneumatic Compression Devices for Venous Thrombosis Prophylaxis in Trauma Patients. Presented at the annual meeting of The American Thoracic Society, May 1994.



SPECIALTY MEDICAL PRODUCTS, INC.

- Billing:** 50 Pennsylvania Ave Malvern, PA 19382
- Shipping:** 187 Pennsylvania Ave Malvern, PA 19382
- Phone:** (888) 698-1522 Fax (610) 644-3992
- Email:** info@smpcares.com
- Website:** www.smpcares.com



3614 Kennedy Road, South Plainfield, New Jersey 07080 • phone: 800-993-9013 • fax: 908-822-1489
Web Site: www.ctcdvt.com

© Copyright Compression Therapy Concepts, all rights reserved. ®Vaso Press is a registered trademark of Compression Therapy Concepts