



SPECIALTY MEDICAL PRODUCTS, INC.

Billing: 50 Pennsylvania Ave Malvern, PA 19382 Shipping: 187 Pennsylvania Ave Malvern, PA 19382

Phone: (888) 698-1522 Fax (610) 644-3992

Email: info@smpcares.com Website: www.smpcares.com

extriCARE® 3600

Negative Pressure Wound Therapy System

Product Overview

The extriCARE® 3600 Negative Pressure Wound Therapy (NPWT) device is a portable, battery powered pump designed for use in hospitals, wound care centers, and other medical facilities. The extriCARE® 3600, used in conjunction with our anatomically fitted line of wound dressings, may promote wound healing through the drainage and removal of wound exudates, infectious material, and tissue debris from the wound bed using continuous and/or intermittent negative pressure.

Features and Benefits:



Versatile: High-capacity, high flow pump, offers maximum care.



Powerful Long lasting rechargeable battery may last for up to 24 hours of treatment.



Easy to Use: Simple operation and user friendly interface.

Device Specifications:

Weight:	2.87 lbs (1.3kg)
Battery Type:	Lithium (rechargeable)
Vacuum Modes:	Continuous or intermittent
Charging Time:	< 5 hours
Dimensions:	Length: 6.7" (17cm)
	Height: 5.1" (13cm)
	Depth: 4.3" (11cm)
Input:	100-240 V
	50/60 Hz
	1.5 A



Item Code: 3600 LB38.0005 Rev B 20160115



extriCARE® Dressings

Negative Pressure Wound Therapy System

Dressings Overview

Our anatomically fitted dressings come in a variety of sizes to fit most needs. They feature a non-adhesive contact layer which can be applied directly to the wound bed. This layer is centered on an adhesive film which creates a tight seal, and may promote healing. A drainage tube connects each dressing to a canister connected to our pump to remove wound material. Our foam dressings come in two sizes in a kit that includes adhesive film, connector tubing, and a suction bell.

Features and Benefits:



Versatile: Available in many sizes, can be used with all extriCARE® pumps.



Simple: Easy to use, quick application



Effective: Creates a tight seal, important for creating negative pressure for treatment of wounds.







Sacral EC2400-S



Small Oval EC2400-SO



Foot EC2400-F



Mini Oval EC2400-MN



Large Foam Dressing EC-FOAM-L



Large Oval EC2400-LO



Small Foam Dressing EC-FOAM-S

LB30.0004 Rev H 20160115