



**SPECIALTY MEDICAL PRODUCTS, INC.**

**Billing:** 50 Pennsylvania Ave Malvern, PA 19382

**Shipping:** 187 Pennsylvania Ave Malvern, PA 19382

**Phone:** (888) 698-1522 Fax (610) 644-3992

**Email:** info@smpcares.com

**Website:** www.smpcares.com

## Model 1850 Knee CPM™



### *Powerhouse of Fundamentals*

The Furniss Corp Model 1850™ Knee CPM combines durability with ease of operation. Each device is made of mild steel construction and built to withstand the rigorous hospital and home health environments.

The full featured CPM offers an easy to comprehend digital display pendant, bariatric accommodations and space-saving tripod position. The Model 1850™ Knee CPM provides proper anatomical therapy, accurate range of motion and ultimate comfort for all patient types. Attention to detail and craftsmanship that goes into constructing each machine by hand in the USA.

### **CPM Devices + Softgoods**

*Manufactured for Providers All Over the World*

US 800.273.5233 Fax 614.871.9353

KneeCPM.com

### Technical Specifications

Knee Flexion: 110 degrees  
Knee Extension: -5 degrees

Speed: 1 to 3 degrees/second  
Pause: 0 to 30 seconds

Limb Length: 20 to 42 inches  
Calf Length: 11 to 24 inches  
Thigh Length: 9 to 18 inches

Power Supply: 110v  
Unit Weight: 25 lbs.  
Ship Weight: 38 lbs.  
Length: 28 inches



# Model 1850 Knee CPM

## Pediatric Athletic Bariatric

Proper anatomical therapy, accurate range of motion and ultimate comfort for all patient types

- Lightweight (25 lbs.) and Portable (28 in.)
- Constructed with mild steel and a dual ball screw assembly—exclusively
- Circuitry located in base of unit instead of hand control for fundamental protection
- Easy to understand hand-held control features a bold LED digital display
- Patient lock out feature prevents inadvertent changes to settings
- Compliance meter records treatment time
- Features safety stop if motion is obstructed
- Universal left/right design
- Quick conversion from adult to pediatric—no additional accessories needed
- Tripod position offers easy storage option and reduces contamination
- Axis to axis motion, skid pad and bed stabilizing system prevent migration
- Knee CPM softgoods are constructed with deep pile synthetic sheepskin and high tensil strength Velcro to ensure patient comfort and proper positioning during therapy. (order #1828-1)
- Two-year extended manufacturer warranty
- HCPCS E0935



### Furniss Corporation Defined

As a US based medical manufacturer, The Furniss Corporation builds Continuous Passive Motion (CPM) devices for the upper and lower extremities. The durable, user friendly devices are built to withstand the rigorous hospital and home care environment. Furniss Corp also manufactures patient softgoods and offers a comprehensive technical service and repair center. Expect same-day shipping, immediate access to educated product professionals and strong product lines with exceptional performance. Visit the website at [KneeCPM.com](http://KneeCPM.com) to experience the interactive tour.

### Indications

Total knee replacement  
ACL reconstruction  
Manipulation under anesthesia  
Ligament reconstruction  
Meniscus repair  
Prosthetic replacement  
Prolonged joint immobilization

Class I Medical Device

CPM Devices + Softgoods  
US 800.273.5233  
[KneeCPM.com](http://KneeCPM.com)



**SPECIALTY MEDICAL PRODUCTS, INC.**

**Billing:** 50 Pennsylvania Ave Malvern, PA 19382  
**Shipping:** 187 Pennsylvania Ave Malvern, PA 19382  
**Phone:** (888) 698-1522      **Fax (610) 644-3992**  
**Email:** info@smplecares.com  
**Website:** www.smplecares.com

## CPM Patient Softgoods



The Furniss Corporation manufactures the custom designed Softgoods for the Phoenix Series Knee CPMs and the DC2480 Knee CPM. We have put the same considerations for quality, fit and function in our softgoods as we have in the design of our Continuous Passive Motion therapy equipment.

We believe that the quality of the softgood is as important as the CPM in delivering patient comfort. We use high quality, deep pile Winning Grey sheepskin and high-tensile strength Velcro brand strapping in our manufacturing. Our Softgoods are designed to maximize the troughing of the leg that is unique to our open carriage CPM design. This allows for more

flexibility in positioning of the patient which increases comfort and the overall quality of the knee therapy.

Furniss Corp ships a free set of softgoods with every Phoenix and DC Model knee CPM. Maintain the quality and re-order the Softgoods made by the company building the CPM.

*Re-Order Information: Part 1828-1 Furniss Softgoods Kit*

### ***Softgoods for the Knee CPM Industry***

Furniss Corp CPM Softgoods are an excellent choice for many competitor CPM devices. The design provides exceptional patient comfort for superior quality and affordability to the healthcare provider.