

Chair-Air™ Series

Alternating Pressure Chair Cushion Systems

CHAIR-AIR™ 9700GR - GERIATRIC RECLINER MATTRESS AND PUMP SYSTEM

This unique alternating pressure mattress system allows for patients who are at high risk for pressure sores to be placed in a variety of mobile devices. Patients are no longer confined to bed for the whole day due to risks of sitting without proper pressure redistribution. Chair-Air™ systems can be integrated into a geriatric chair, recliner or standard high-back wheelchair. Twenty bladders gently alternate, providing low interface pressures in addition to comfortably supporting the patient. The Chair-Air™ 9700GR system includes the 9701 alternating pump. Chair-Air™ meets California Technical Bulletin #117 for fire retardancy. An optional fire barrier is available to meet the most stringent fire safety requirements, when requested.



Recliner not included

CHAIR-AIR™ 9700C - ALTERNATING AIR WHEELCHAIR CUSHION AND PUMP SYSTEM

Six individual bladders alternate to provide the patient a superior pressure redistribution seating surface. The zippered non-skid, stretch-knit cover is anti-microbial and fire-retardant. The Chair-Air™ 9700C includes the Chair-Air™ 9701 alternating pump.



Shown in cover

CHAIR-AIR™ 9700CAF - ALTERNATING AIR WHEELCHAIR CUSHION AND PUMP SYSTEM

Four individual bladders alternate to provide the patient a superior pressure redistribution seating surface while the unique molded design, with slight abductor pommel, stabilizes the hips and keeps the patient in proper mid-line alignment. The zippered non-skid, stretch-knit cover is anti-microbial and fire-retardant. The Chair-Air™ 9700CAF includes the Chair-Air™ 9701 alternating pump.



Shown in cover

CHAIR-AIR™ AC/DC 9701 - ALTERNATING PUMP

The portable Chair-Air™ AC/DC Pump allows your patient to be totally mobile. A rechargeable battery backup keeps this system working up to 6 hours per charge, based on the patient's weight.



For use with Chair-Air™, Geriatric Recliner Mattress & Chair-Air™ Wheelchair Cushion



"Setting the Standard"

CHAIR-AIR™

SPECIFICATIONS

CHAIR-AIR™ GERIATRIC RECLINER MATTRESS - MODEL 9702GR

Dimensions: 19" x 68" x 2.5 - 3.5"
Weight: 6 lbs.
Air Cells: Multiple - 2" - 3.5" sacs

Recliner not included

Mattress Cover:

Low shear, zippered multi-stretch Vyvex™ III, moisture / vapor permeable top cover, non-skid bottom. Cover is anti-microbial and meets California Technical Bulletin #117 for fire retardancy.




CHAIR-AIR™ WHEELCHAIR CUSHION - MODEL 9702C

Dimensions: 18" x 16" x 3.25"
Weight: 4 lbs.
Air Cells: 6 - 3.25" sacs



CHAIR-AIR™ WHEELCHAIR CUSHION - MODEL 9702CAF

Dimensions: 18" x 16" x 3"
Weight: 3 lbs.
Air Cells: 4 - 2.5" sacs



WHEELCHAIR CUSHION COVER:

Zippered low shear stretch-knit with non-skid bottom. Cover is anti-microbial and meets California Technical Bulletin #117 for fire retardancy.

PUMP FOR USE WITH GERIATRIC RECLINER & STANDARD WHEELCHAIR CUSHION

Dimensions: 14.5L" x 5W" x 12.5H"
Weight: 16 lbs.
Electrical: A/C 117V/60 hz D/C 12V 7 amp/hr
Listings: UL CE
Cycle Time: 10, 15 or 20 Minutes
Air Flow: 4 - 6 Liters
Ground Resistance: 0.1 ohms maximum



CHAIR-AIR™ AC/DC 9701 ALTERNATING PUMP

WARRANTY

Pumps: 1 year limited
Cushion: 1 year limited
Cover: 6 months limited



7-11 Suffern Place, Suffern NY 10901
 Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115
 www.bluechipmedical.com

"Setting the Standard"

©Blue Chip Medical Products, Inc. 2015

- Alternating Pressure Mattresses
- Alternating Pressure Pads & Pumps
- Low Air-Loss Mattress Systems
- Lateral Rotation Mattress Systems
- Gel Mattresses
- Gel Mattress Overlays
- Static-Air Mattress Overlays
- Water Mattress Overlays
- Wheelchair Cushions
- Patient Lift Slings
- Side Rail Bumpers
- Cart Covers
- OR Replacement Table Pads
- Replacement Stretcher Pads
- Cubicle Curtains
- Hospital Pillows